Naipo

User Manual

Wireless Eye Massager with Heating & Music

Model: MGE-1S3

Contents

English

Deutsch

Français

Español

Italiano

01-03

04-06

07-09

10-12

13-15

Production Device

The eye massager includes main body, USB cable and user manual.

Operation Guide

Preparation before use

① Full Charge

② Our Company's Power Adapter

Operation

On: Press the power button momentarily, it works after a beep.

Off: Press the power button for two seconds, it stops after three beeps.

Mode Operation

Press the power button momentarily to turn on.

Working Situation:

The appliance is in sleep mode by default.

Press the power button momentarily to change modes. It shows MED/HARD/SLEEP in turn.

Press the power button steadily to turn off.

SLEEP MODE: 10 minutes

Air compression: on

Heat: on

Vibration: off

MED MODE: 15 minutes

Air compression: off

Heat: on

Vibration: on

HARD MODE: 20 minutes

Air compression: on

Heat: on

Vibration: on

Product Functional

Point Massage

Air Pressure

Heating

Soothing Music

Vibration

Product Functional Configuration

Power Button

Music Button

USB Interface / Charging Interface

Basic Function of Production

Air Pressure	Heating	Time	Vibration	
Mild / Intense	40 °C~45 °C	Stream Wave Sound Clocks Airplop	10 minutes 15 minutes 20 minutes	Default Mode Unchangeable

Operation of Music

- Press the music button momentarily to turn on.
- Press the music button momentarily to change songs.
- Press the music button steadily to turn off.

Notes and Statements

① When battery is low, please charge it with external power device. The massager can't be used when charging. ( Battery symbol shows full means workable, while it flashing means low battery. )

② On different modes, eye massager will show corresponding modes and related functions.

③ All modes are fixed, time and function can't be adjusted.

④ The machine will automatic shutdown in low battery, please charge it with external power device to make it full.

⑤ The machine will beep while adjusting modes.

Cautions

Please read this user manual carefully before using this product. Naipo is not responsible for any results of improper usage or violation of these instructions.

① Please do not use the eye massager when driving.

② Please do not use the eye massager in excessive high temperature or moisture.

③ Please do not use the eye massager when you have eye diseases such as glaucoma, cataract, retinal detachment and so on.

④ People who have an eye surgery should consult a medical professional before use. Stop the use immediately if you have any discomfort.

⑤ Patients who suffer from severe stroke, high blood pressure and low blood pressure should consult a medical professional before use. Please stop the use immediately if you have any discomfort.

⑥ People who suffer from severe slow-moving, physical and mental discomfort and mental retardation should use with caution under the supervision of attendants. Please stop the use immediately if you have any discomfort.

⑦ Please do not use eye massager after drinking.

⑧ Please ensure no contact lenses before use.

⑨ Please turn it on after wearing eye massager, so as to avoid damages of eyebag over-pressure.

⑩ Please use adapter of this massager and ensure the safety of electricity use.

⑪ This product is not suitable for diagnosis, treatment and other commercial uses.

⑫ Please wipe it with disinfection cotton soaked in neutral cleaner instead of water.

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Shenzhen Pirlong E-Commerce Co., Ltd.