

印刷：单黑双面印刷


材质：105g铜版纸

尺寸：70*120mm

折页方式：风琴折



User Manual



Rechargeable TENS Unit

Model: MGEP-AS1080

Introduction

Thank you for purchasing the Naipo Rechargeable TENS Unit. Please read this user manual carefully and keep it in a safe place for future reference. If you need any assistance, please contact our support team with your product model number and Amazon order number.

Included


- 1 x TENS unit
- 3 (pairs) x Large pads
- 2 (pairs) x Small pads
- 1 x Storage plate
- 2 x Connecting wire
- 1 x USB cable
- 1 x Charger
- 1 x User manual
- 1 x Storage pouch

Product Specifications


Model	MGEP-AS1080
Name	Naipo Rechargeable TENS Unit
Output frequency	0-452HZ
Output pulse width	245US
Output voltage 0.5K load	0-50
Output pulse circuit	2
Built-in battery	Lithium battery (180MA)
Output voltage of the adapter	5V
Input voltage	AC110-220
Working time	0-60 minutes
Environment temperature	-20°C-40°C
Relative humidity	≤85%

- 01 -


Product Diagram



- 01. Two output ports
- 02. USB port
- 03. LCD screen
- 04. Battery indicator
- 05. ON/OFF switch
- 06. Intensity adjustment buttons
- 07. Mode selection button
- 08. Time adjustment button
- 09. Lock button / Secondary mode selection button
- 10. Storage plate
- 11. Charger
- 12. USB cable
- 13. Massage pad
- 14. Connecting wire



- 02 -



- ON/OFF Switch: Slide the switch to the "ON" position to turn on the product with the screen displaying all things, and the product will shift to standby mode after 2 seconds. The default working time is 20 minutes (when time is up, the product will shut off automatically).
- The product is equipped with 6 primary modes (16 secondary modes for primary modes): Knocking (with 4 secondary modes), Acupuncture (with 3 secondary modes), Cupping (with 3 secondary modes), Tapping (with 2 secondary modes), Shiatsu (with 2 secondary modes), Massage for Sensitive Parts (with 2 secondary modes).
- The product is equipped with 5 function buttons (intensity decreasing button / secondary mode selection button / mode selection button / time adjustment button / intensity increasing button):
 - Intensity decreasing button: Decrease an intensity level with each press of this button (weakest intensity: level 1) Secondary mode selection button (the round button): Press to choose different secondary modes; Long press for 3 seconds to lock the function buttons and long press for 3 seconds again to unlock.
 - Mode selection button: Press to interchange among 6 primary modes.
 - Time adjustment button: There are 6 timer settings: 10/20/30/40/50/60 minutes.

- 03 -

Shenzhen Pinlong E-Commerce Co., Ltd.

Web: www.naipocare.com

E-mail: support@naipocare.com

Add: Huanan City Electronic Trading Building P07

Room 120, Longgang District, Shenzhen, China.



- Intensity increasing button: Increase an intensity level with each press of this button (highest intensity: level 20)

Product Operation Guide

- Install the pads on the connecting wire, then remove the protective film.
Note: Adding a few drops of water on the pads will help prolong their service life.
- Apply the pads on body areas where massage is needed. Connect the connecting wire to one of the output ports on the product. For additional relief, repeat the procedure with another port being used to apply another 2 pads.
Note: Make sure the pads on the skin are not overlapping and that your skin is free of any dirt, oil, or lotion.
- Slide the ON/OFF switch on the top to the "ON" position to turn on the product.
- Press the mode selection button to choose among primary modes, and press the central round button to choose among secondary-level modes. Remember that the intensity would shift to the lowest level every time a mode is changed.
Note: The central round button can also be used as a lock button (Long press for 3 seconds to lock the function buttons and long press for 3 seconds again to unlock)
- After selecting a mode, increase the intensity via the "+" button and decrease the intensity via the "-" button.
- The recommended time to use the product is 10-60 minutes. Increase 10 minutes with each press of the time adjustment button. When time is up, the product will shut off automatically.

Important Safeguards

- There should be a good connection between the pads and the connecting wire.

- 04 -

- The pads should be applied on the proper areas on the skin.
- Connect the connecting wire to the main body via the DC hole.
- After checking the connection, turn on the product; Please turn off the product before removing the pads from the skin.
- Do not stick a pad to another to avoid short circuit.
- Do not use the product on or near the heart, the head, the neck, the privates, the injured areas, the chest, or the mouth.
- Do not use the product while driving, operating machinery, or undertaking exercise with unconscious muscle shrinkage to avoid possible dangers or injuries.
- People listed below should not use the product:
 - People who have suffered from an acute disease or trauma, or people who have experienced a surgery in the past 6 months.
 - People with cardiac dysfunction.
 - People with epilepsy.
 - People with high blood pressure.
 - People with a high fever.
 - People with cancer.
 - People with an infectious disease.
 - People with a malignant tumor.
 - People who have experienced a surgery recently.
 - People with skin problems.
 - Pregnant women.
 - People with undiagnosed pain syndromes.
- Do not touch the pads when the product is being turned on.
- Do not use the product if you are equipped with a pacemaker or other medical medical apparatus like an artificial lung or a respirator.
- Keep the product beyond reach of children.
- Stop using the product immediately is pain if felt.

- 05 -

Questions & Answers

Problem	Cause	Solution
The product does not work.	There exists a loose connection or improper operation.	Please check whether the product has been connected properly. 2 pads should be used simultaneously and the protective film should be removed before usage.
The pads could not be attached to the skin well.	There exists dust or dirt on the pads, which might weaken the electrical conductivity.	Please keep the skin clean without sweat before each use. Attach the protective film back on the pads promptly after usage. Do not use under direct sunlight or in humid place.
The intensity from one pad is stronger than that from the other.	It's a normal phenomenon, for different areas of your body will react differently.	Nothing needs to be done. Make sure that the pads are moist and that they contact the skin well.

- 06 -

During the massage, a painful burning sensation is felt on the skin or the stimulation becomes weakened.	The pads are not attached firmly on the skin.	Apply a few drops of water to the adhesive surface of each pad and make sure that the pads are attached firmly on the skin during the use.
	The pads are too dry.	
The product is on, but no sensation is felt from the pads.	Both pads are not firmly attached on the skin.	Make sure both pads are firmly attached on the skin.
	There is a loose connection.	Make sure all connections are in good condition (the connection between the main body and the connecting wire; the connection between the connecting wire and the pads).
	The intensity level needs to be increased.	Increase the intensity level.
The pads could not be stuck on the skin even after cleaning and moistening them.	The pads need to be replaced.	Replace with new pads.
The product could not be turned on or the screen is dim.	The battery is running low.	Recharge the battery.

- 07 -

Product Maintenance

- After using the product, turn off it, and remove the pads. Place the product in a cool and dry place for next use.
- After usage, please attach the protective film back on the pads and store the product in the storage pouch; Adding a few drops of water on the pads before and after usage will help prolong service life.
- If the product is not used regularly, charge it at least in every 3 months; Charge the product by connecting it to a charger or a computer via the USB port and the USB cable.

Recommend Massage Areas



- 08 -

Warranty & Customer Support

For questions, assistance, support, or warranty claims, please contact us via the email or the address below.

Web: www.naipocare.com

E-mail: support@naipocare.com

Address: Huanan City Electronic Trading Building P07 Room 120, Longgang District, Shenzhen, China.

*Note: Products purchased directly from Naipo benefit from our after-sales service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.

- 09 -